



VIVID PS Developer

Positive Plate Developer for PS plates

Product Usage

Developer Concentrate for Positive Analoge Plates - ready to use

universal for all popular positive working PS plates

Characteristics

- Highly concentrated (1 part PS Developer + 5 parts water) - cost efficient
- Good buffer capacity in the plate processor - long lasting
- Low replenishing rate with Replenishers - cost efficient.
- No or very low residues in the processor tank
- Universal in use - for all popular positive working plates

Application

Dilute 1 part PS Developer Concentrate with 5 parts of water and mix well.

processor adjustments:

Temperature: 22°C-26°C, start with 24°C

Processor speed: 20-40 seconds, start with 30 seconds

Replenishing rate: 15-40ml/m², start with 20ml/m²

Packaging

5/20 liter

Vivid Manufacturing Co. Pvt. Ltd., 14/7, Site-IV, Industrial Area, Sahibabad-201010,U.P.
E-mail : info@vividchem.com WebSite: www.vividchem.com